

# Product datasheet

## Gymna Mobile

*Practical and nice looking trolley for ambulant use of your therapy devices. Matches the Gymna 200 series perfectly.*

### Standard features

- Top shelf in fixed position, 90 cm high.
- Lower shelf is adjustable from 50 - 80 cm in six steps.
- Handy drawer underneath the lower shelf.
- Integrated bottle holder in the top shelf, very convenient for contact gel or massage lotion. The bottle can also be stored upside down.
- Detachable cable and/or vacuum cup holder for tidy storage of your accessories.
- Sturdy central column that allows hiding the power cables.
- Equipped with four swivel wheels for easy manoeuvrability.
- Two front wheels can be locked to hold the trolley in position.

The Gymna Mobile is delivered in a Do-It-Yourself package and requires little time for set-up.

### Ordering data

320 804 - Gymna Mobile equipment trolley

### Specifications

Height top shelf:	90 cm
Total height:	94 cm
Width x depth:	56 x 53 cm
Height 2 <sup>nd</sup> shelf:	50 - 80 cm
Weight:	± 20 kg
Max. load per shelf:	15 kg
Colour:	white RAL 9016
Shelves:	laminated (HPL) plywood
Packaging:	83 x 58 x 21 cm, 23 kg

